EN-ECG Cable& lead wire

Direction for use

Operation principle:

ECG is a presentation of physiological activity .A sequence of physiological change would generate during the heartbeat. These physiological change will transit to the human-body surface. The ECG lead wire electrodes which are connected with the machine can catch the electrical signal. The amplified and processed signal will be converted to waveform and digital signal which will be displayed on the patient monitor.

Product description:

Product name: ECG Cable&lead wire

Product model:

2307P,2307S-I,AA-2375,AA-2385,LL5-90P,SL5-90P,AAB5-90DS,MQ5-90DS,E10R-HL,2308P,2308S-I,AA-2575,AA-2585,SM5-90P,MN5-90P-I,NK-90DS,E10-FD2-B,2309P,2309P-I,D-1574-I,D-1307,AA3-90P,SM3-90P-I,E10R-HP-B/I,2310S,2311S,D-1308,D-1309,CL3-90P,AA5-90S,E10R-MQ-B/I,2312P,2320S,D-1310,D-1311,D3-90P,CL5-90S,E10-NK2-N,2330P,2332P,DX-2355,DX-2395,D5-90S,DX5-90P,E10R-SH2-P,2335S,2340P,DX-2595,SL-2540,NECA3-90P,D7-90S,E10R-FD1-B,2340S,2341P-I,HP-2385,HP-2585-I,HP5-90P,MN3-90P,2355P,2374P,LL-2340,MQ10-LB,HP3-90S,MQ5-90P,2375S,2376P-I,MQ-2586,LL-2540,D10-90S,2385P,2386S,NK-2507,NK-2509,2393P,2395P,NK-2540,D-1555,2396P-I,2397P,D-1540,2398P,2507P,2508P,2510S,2511P,2512S,2530P-I,2532P,2540P,2540S,2541P,2555P,2574P,2575P,2576S,2585P,2586P-I,2595P,2596P,2597P,2598S,2599S,HPA5-90DS,AT3-90DS,E10-BK1-B,E10-FD1-B,E10-KZ-B,E10-NK1-N,E10R-SH1-B,E10-MT-B.

Product specification: one piece fixed ECG cable、ECG trunk cable、lead wires

Primary structure: It composed by jacket plug/cable/YOKE/the plug of YOKE/electrode connector etc.

Application:

It's suitable for multi-parameter monitor/electrocardiogram machine/exercise treadmill or defibrillator ,transit the bio-electricity signal of human-body.It's suitable for all population and common place.

Environmental requirements:

Environment temperature: 0°C ~40°C

Relative humidity: 15% ~85%

Atmospheric pressure: 86kpa ~106kpa

Package and storage:

Every ECG cable should be packaged singly, it should be stored in drafty indoor without corrosive gas.

The storage conditions are below:

1. Environment temperature: -10°C ~+40°C
2. Relative humidity: ≤80%

Safety:

Meet with the regulation of <medical electrical equipment first part: general safety requirements>.

Installation and use:

Connect the ECG lead wire with the adaptive patient monitor, and check whether the operation is the same as the operation instructions of operation manual.

The ECG electrode of lead wire should be connected on the patient in order, the lead wire should be connected separately to avoid entwine, otherwise it may cause interference during the monitoring, such as heart rate fluctuation, waveform fluctuation etc.。

The ECG electrode of lead wire can be divided into the following kinds:


- 1), grabber end
- 2), snap end
- 3), little grabber end

- 4), 4mm banana end
- 5), 3mm needle end

The connection type of electrode on the human-body as below:

Three leads AHA

Five leads AHA

Ten leads AHA

AHA VS IEC comparison table:

AHA function	IEC function	human body parts
White RA	red R	right arm
Black LA	yellow L	left arm
Red LL	green F	left leg
Green RL	black N	right leg
Brown V	white C	single chest lead
Red V1	red C1	forth intercostal of breastbone right edge
Yellow V2	yellow C2	forth intercostal of breastbone left edge
Green V3	green C3	the middle point of the line between C2 and C4
Blue V4	brown C4	the crosspoint of midclavicular line and fifth intercostal
Orange V5	black C5	the crosspoint of the same level of left anterior axillary line and V4
Purple V6	purple C6	the crosspoint of the same level of left midaxillary line and V4

Cleaning:

The cable should not be connected with the patient monitor when it is be cleaning or disinfecting.

It should be cleaned or disinfected before used with another patient.

It could be cleaned the surface of ECG electrode with 70% isopropyl alcohol or ethyl alcohol solution.

Use a clean/dry sponge cushion dip into the clean solution in advance when cleaning. Scrub all the surface of the electrode and the cable with this cushion. Use another clean/dry sponge cushion dip into the sterile or distilled water in advance. Scrub all the surface of the electrode and the cable with this cushion. Finally dry all of the surface of the electrode and the cable with a clean/dry sponge cushion.

Caution : the pins of jacket plug should not contacted with the clean solution, otherwise it will cause damage to the patient monitor and the human body.

Warranty period/period of validity:

Warranty period:6 month

period of validity: N/A

Replacement and discard:

Replacement: the breakage of cable/the damage of electrode/the damage of outer jacket cable/ expiration of period of validity of such adverse circumstance, please replace the new ECG cable in time.
Discard: The discarded ECG lead wire is the medical waste material, please seek the professional department for recycling.

Contraindication/notice/caution/indicative:

Contraindication: N/A

All the lead wires which are connected with the patient, the incorrect connection or failure may cause reading error. High temperature or electromagnetic wave interference and the dangerous voltage source may cause inaccurate reading also.

The lead wire should be checked that whether the connection with the equipment is reliable every 2 hours.

The cable should be placed correctly to avoid the cable entwine the patient cause suffocate.

The lead wires should be installed and placed correctly according to the requirements of use manual of equipment.

It strictly prohibited that it should not be used with the MRI equipment at the same time, otherwise it

may cause patient burn.

When it is connected with the defibrillator. The electrodes on the patient should keep distant as far as possible with the electrode of defibrillator, otherwise it may cause burn between the electrode and the contact segment of patient when the electric discharging.

The stretch or bent to the wire should not be overexertion, otherwise it may cause the damage of internal cable and the lost of the shield properties or the inaccurate reading.

The metal portion of the lead wire would keep dry, otherwise it may cause failure.

Refer to the operation manual of patient monitor please.

Symbol explanation:

Caution	Series number	Manufacturer	Non-sterilization	Latex free
EU representative	Temperature limitation	Humidity limitation	See the instructions	Don't throw Away-Reusable
Expiry date	BF type equipment	CF type equipment	Dispositifmedical	CH representative

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