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**Leadwires**

**Verbindungskabel**

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EN

## Instruction of use and storage for leadwires

The instruction is concerned with leadwires type: PD, W.KAB.PD, PG, W.KAB.PG

### 1. General description

Leadwires are intended to transmit electrical signals received from human organism to devices processing these signals - ECG apparatus, unassisted or with consolidated cable.

Leadwires are non-sterile class I medical products.

### 2. Leadwires characteristics

Type of applied wire, its length, type of endings on patient's side, usage of resistance or inductive elements are adjusted to specific type/model of apparatus or consolidated cable with which the leadwire is to cooperate.

Marking of the leadwire, method of its connection, purpose of individual tips and outlets, and other data essential for correct conduction of examinations, is specified in the instruction of use of the device to cooperate with which the leadwire is intended.

Leadwires are intended for use with undamaged patient's skin.

Leadwires should be used exclusively by qualified medical personnel.

### 3. Principles of use

When using leadwires, the following principles should be observed:

- leadwires should not be excessively bent (bend radius should not be lower than 60 mm),
- do not pull or jerk the leadwire,
- place the leadwire so that there would be no risk of mechanical damage, e.g. by stepping on, running over, placing a heavy/sharp object on it, cutting etc.,
- hold the connector cover (socket, pin, clasp, latch ) when connecting and disconnecting the leadwire,
- protect the leadwire connectors against suffusion.

### 4. Installation

Connect the cable pin to the appropriate consolidated cable/ECG apparatus, turn on the device and verify the correctness of functioning according to monitor's instruction of use.

Pay attention to proper choice of the leadwire for appropriate consolidated cable/ECG apparatus.

## 5. Recommended cleaning, disinfection

These activities should be realized every time before use of the device and leadwire before the examination of the next patient. Before starting the cleaning or disinfection the leadwire should be disconnected from the device with which it cooperates.

### Cleaning

Wash the leadwire by gently wiping it with soft cloth moistened with warm water with soap or mild detergent. The leadwire should be left to dry before use.

Do not use corrosive substances or other chemical preparations.

When cleaning, do not soak or immerse the leadwires, do not wash under running water, and do not use excessive strength which may rupture the internal conductors, leading to leadwire damage. Continuous bending of the leadwires during use and cleaning may rupture the internal conductors. Damage of this type is not covered by the guarantee.

Avoid contact with strong solvents: aromatic, chlorinated, ketone, ether and ester. Do not allow for contact of cleaners with tips from the electrode side and pins from cable / ECG apparatus side.

### Disinfection

Disinfect leadwires by means of preparations dedicated to disinfect medical devices made of plastics, based on active ingredients: isopropyl alcohol (70% concentration), glutaraldehyde, amines and quaternary ammonium compounds. The disinfection should be carried out in accordance with the instructions of the manufacturer of the preparation, and then leave the leadwire to dry. Do not allow for contact of preparat with tips from the electrode side, and pins from cable / ECG apparatus side.

## 6. Packaging and storing

Leadwires are packed individually. In order to maximize essential characteristics of the products, they should be stored in the original packaging, in accordance with provided storing conditions.

Storage conditions:

- temperature: -25°C to +55°C,
- relative humidity: 10% to 85%, do not allow for steam condensation
- protect against sunlight, moisture and dirt.

Leadwires in use should be stored in a dry place.

When not used, leadwires should be loosely coiled, without sharp bends.

## 7. Environment conditions during work

- temperature: +5°C to +45°C
- relative humidity: 30% to 60%
- atmospheric pressure: 700 to 1060 hPa

## 8. Warnings

Not observing the principles provided above may break electric circuit or lead to short-circuit. These signs are easy to identify by the ECG apparatus service staff. Damaged leadwire should be disconnected and replaced with an operational one. Leadwire should be used with fully functional devices, in the event of damage to the leadwire or ECG apparatus, stop the test immediately. Use only in rooms equipped with an efficient electrical and grounding installation.

Leadwires are intended for use with specified types of ECG apparatuses. The user is responsible for verification, before its use, of leadwire compatibility with the device. Incompatibility of these elements may lead to reduced accuracy and efficiency of measurements (or lack of them), and in extreme case - to damaging of the product. In order to ensure their compatibility, check the equipment cooperating with a given device in its instruction of use.

It is not recommended to use leadwire in excitable patients because movement negatively influences realized measurements. Relocation and positioning of the patient may cause disconnection of the leadwire from ECG electrodes and lack of readings.

Strong sources of electromagnetic field, such as those emitted by electrosurgical equipment, may negatively influence the correct functioning of the leadwire. Do not use the leadwire during magnetic resonance imaging (MRI).

Allergic reactions at contact of leadwire with patient's undamaged skin are not observed. However, such reactions may not be excluded for persons particularly prone to allergies, the manufacturer are not known for this type of event. Improper cleaning or disinfection (or lack of them), and disregard for hygienic procedures may result in transfer of micro-organisms and cause infections.

May used during defibrillation, leadwires should be used in conjunction with cable / ECG equipment with defibrillation protection.

Do not use leadwire suspected of having mechanical and/or electrical defects, or packaging of which has been damaged. Suffused leadwire should not be used.

Do not change or modify the leadwire in any way. Such actions may influence correctness of functioning and accuracy of received readings.

Utilisation of used leadwire should be done in accordance with binding provisions.

## 9. Guarantee

A life time of the product is 10 years from the date of manufacture.

The guarantee period covers 12 months from the date of sale.

The general terms and conditions of the guarantee are described in a separate document available at the customer's request and at [www.sorimex.eu](http://www.sorimex.eu)

The guarantee does not include defects resulting from improper use and storage.

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