1. Product Introduction

The ECG trunk cable serves as a crucial connection between the ECG leadwires and the main ECG device. It aggregates the signals from multiple leadwires and transmits them to the device for processing, playing a vital role in the ECG monitoring and recording system. Our ECG trunk cables are engineered for stable signal transmission and compatibility with a range of ECG setups.

2. Product Components

- **Trunk Cable:** The main cable is constructed from high - quality materials to ensure efficient signal transmission. It has a gray outer insulation for identification and protection.
- **Connectors:** One end of the trunk cable is designed to connect to the ECG leadwires (with multiple ports to accommodate different leadwire connectors). The other end connects to the ECG device. The connectors are designed for secure and reliable connections, with proper pin configurations to match the device and leadwires.
- **Labels and Markings:** Each trunk cable is labeled with a model number, serial number, and relevant safety and compliance markings.

3. Installation Process

3.1 Preparation

- Turn off the ECG device and any connected leadwires before installing the trunk cable.
- Inspect the trunk cable for any visible damage, including cuts, frayed cables, or broken connectors. Do not install a damaged cable.

3.2 Connecting to ECG Leadwires

- Identify the appropriate ports on the trunk cable for connecting the ECG leadwires. The ports are usually designed to match the connectors of the leadwires (e.g., if the leadwires have a 5 - pin connector, insert it into the corresponding 5 - port on the trunk cable).
- Gently insert each leadwire connector into the trunk cable port until it is fully seated and secure.

3.3 Connecting to ECG Device

Locate the ECG device port that is compatible with the trunk cable connector. Align the connector and insert it firmly into the device port. Ensure a proper connection to avoid signal loss.

3.4 System Activation

Power on the ECG device and the connected leadwires (if applicable). The trunk cable will now transmit the aggregated signals from the leadwires to the device.

4. Maintenance and Care

- **Routine Checks:** Regularly examine the trunk cable for any signs of damage. Check the connectors for bent pins, corrosion, or loose connections. Inspect the cable for kinks, cuts, or wear.
- **Cleaning:** Wipe the trunk cable and connectors with a soft, dry cloth. For disinfection, use a mild medical - grade disinfectant wipe. Avoid using harsh chemicals that could deteriorate the cable insulation or connectors. Do not submerge the trunk cable in water or cleaning solutions.
- **Storage Conditions:** Store the trunk cable in a clean, dry area. Keep it away from extreme temperatures (ideally between 10°C - 40°C) and high humidity. Coil the cable loosely to prevent permanent kinks.

5. Troubleshooting

- **Intermittent Signal or Signal Loss:**

- Check all connections (between leadwires and trunk cable, and trunk cable and device). Reconnect if necessary to ensure a secure fit.
- Inspect the trunk cable for any physical damage that could be causing intermittent contact.
- Test the trunk cable with a different set of leadwires or a different ECG device (if available) to isolate the problem.





- **Incompatibility Issues:**

- Ensure that the trunk cable is compatible with both the ECG leadwires and the ECG device. Check the product specifications and compatibility charts.
- If a new device or leadwires have been introduced, verify that they are within the supported range of the trunk cable.

6. Safety Guidelines

- Never use a damaged trunk cable, as it may compromise signal integrity and pose a risk of electrical hazards.
- Follow the manufacturer's instructions for the ECG device and leadwires when using the trunk cable. Incompatible combinations may lead to inaccurate results or device malfunction.
- Keep the trunk cable away from sources of electrical interference, such as large motors, transformers, or other high - frequency equipment, to ensure accurate signal transmission.
- Adhere to all relevant medical device safety regulations and standards in your region when installing and using the ECG trunk cable.

7. Symbol explanation:

Symbols	definition	Symbols	Definition
	Note investigation, please consult the documentation that came with		Product complies with the EU Medical Devices Directive 93/42 / EEC of identity
	Waste Electrical and Electronic Equipment	Rx ONLY	Federal (U.S.) Law restricts this device to sale by or on the order of a physician
	Storage temperature limits		serial number
	Storage humidity limits	REF	Catalogue number
	Large oxygen pressure limits		Company logo and Do not discard
	Latex-free		EU authorized agent
	Non-sterile packaging		manufacturer
	Please refer to the instruction manual		Manufacture date
		IPX2	Waterproof level

8. Manufacturer and EC Representative:

	Name: Hunan Greatmade Medical Tech Limited Add : NO.101, 1st floor, A2 building, 1st Phase, Changsha E center, No.18 Xiangtai road, Liuyang Economic Development zone, Changsha, Hunan Province, China SRN: CN-MF-000042754 Tel :86-731-83167851 Website: www.gmmmedical.net Email: sales@greatmade.com.cn
	Name : SUNGO Europe B.V. Add: Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands SRN: NL-AR-000000247

9. Revision information

Rev.	Revision description	Revised date
A.0	First edition	2025/8/12
A.1	Includes English, German language	2025/8/14
A.2	Includes English, German, French, Spanish, Portuguese language	2025/8/15

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