## 1. Product Overview

ECG leadwires are essential components for electrocardiogram (ECG) systems. They are designed to transmit electrical signals from the patient's body to the ECG device accurately, enabling the monitoring and recording of heart activity. Our ECG leadwires come in various configurations to meet different clinical needs, ensuring reliable signal acquisition.

## 2. Components and Identification

- **Leadwire Cables:** Made of high - quality, flexible materials to ensure patient comfort and durable signal transmission. The outer insulation is gray, providing good protection and identification.
- **Connectors:** Different types of connectors are available to match specific ECG devices. For example, some connectors have multiple pins for connection to the ECG trunk cable or the device directly. Color - coded clips or electrodes at the patient - end help in correct placement on the body (e.g., red for right arm, yellow for left arm, green for right leg, black for left leg in some standard configurations).
- **Labels:** Each leadwire is labeled with a model number (e.g., MC149 - 51G, MC149 - 31G as shown in the pictures) and a serial number for traceability and identification.

## 3. Installation and Use

### 3.1 Preparation

Ensure the ECG device is powered off or in a standby state before connecting the leadwires.

Check the leadwires for any signs of damage, such as frayed cables or broken connectors. If damaged, do not use.

### 3.2 Connection to ECG Device

Align the connector of the ECG leadwire with the corresponding port on the ECG device or the ECG trunk cable. Gently insert and ensure a secure connection. The connector should fit snugly without excessive force.

### 3.3 Patient Attachment

Identify the appropriate placement of the patient - end electrodes or clips according to the clinical protocol. For example, for standard 12 - lead ECG, attach the color - coded electrodes to the respective limbs and chest positions.

- AHA VS IEC comparison table

AHA function	IEC function	human body parts
White RA	Red R	Right arm
Black LA	Yellow L	Left arm
Red LL	Green F	Left leg
Green RL	Black N	Right leg

Brown V	White C	Single chest lead
Red V1	Red C1	Fourth intercostal of breastbone right edge
Yellow V2	Yellow C2	Fourth intercostal of breastbone left edge
Green V3	Green C3	The middle point of the line between C2 and C4
Blue V4	Brown C4	The crosspoint of midclavicular line and fifth intercostal
Orange V5	Black C5	The crosspoint of the same level of left anterior axillary line and V4
Purple V6	Purple C6	The crosspoint of the same level of left midaxillary line and V4

Ensure good skin contact by cleaning the patient's skin (removing oils, lotions, or sweat) and, if necessary, using electrode gel (not included with the leadwires but may be required for some applications).

### 3.4 Operation

Power on the ECG device and start the monitoring or recording process. The leadwires will transmit the electrical signals from the patient's heart to the device for analysis.

## 4. Maintenance and Cleaning

**Regular Inspection:** Periodically check the leadwires for wear and tear. Inspect the cables for any cuts, kinks, or discoloration. Check the connectors for bent pins or corrosion.

**Cleaning:** Wipe the leadwires and connectors with a soft, dry cloth regularly. For more thorough cleaning, use a mild disinfectant wipe. Avoid using abrasive cleaners or solvents that may damage the insulation or connectors. Do not immerse the leadwires in water or cleaning solutions.

**Storage:** Store the leadwires in a clean, dry, and dust - free environment. Coil them neatly to avoid kinking the cables. Do not store them in extreme temperatures (below 0°C or above 50°C) or high - humidity environments.

## 5. Troubleshooting


- **Poor Signal Quality:**
  - a. Check the patient's skin contact. Ensure the electrodes are properly attached and the skin is clean.
  - b. Verify the connection between the leadwire and the ECG device or trunk cable. Reinsert the connector if necessary.
  - c. Look for any damage to the leadwire cable. If damaged, replace the leadwire.
- **No Signal Transmission:**
  - a. Check if the leadwire is properly connected at both ends.
  - b. Ensure the ECG device is powered on and functioning correctly.
  - c. Test with a known - good leadwire to determine if the problem is with the leadwire or the device.

## 6. Safety Precautions

- Do not use damaged leadwires, as they may pose a risk of electrical shock or inaccurate signal transmission.
- Only use the leadwires with compatible ECG devices. Using incompatible devices may result in malfunction or inaccurate readings.
- Keep the leadwires away from sources of heat, open flames, or sharp objects that could damage them.
- Follow all local medical device regulations and safety standards when using the ECG leadwires.

## 7. Symbol explanation:

Symbols	definition	Symbols	Definition
	Note investigation, please consult the documentation that came with		Product complies with the EU Medical Devices Directive 93/42 / EEC of identity

	Waste Electrical and Electronic Equipment	Rx ONLY	Federal (U.S.) Law restricts this device to sale by or on the order of a physician
	Storage temperature limits	SN	serial number
	Storage humidity limits	REF	Catalogue number
	Large oxygen pressure limits		Company logo and Do not discard
	Latex-free	EC REP	EU authorized agent
	Non-sterile packaging		manufacturer
	Please refer to the instruction manual		Manufacture date
		<b>IPX2</b>	Waterproof level

## 8. Manufacturer and EC Representative:

	<p>Name: Hunan Greatmade Medical Tech Limited  Add : NO.101, 1st floor, A2 building, 1st Phase, Changsha E center, No.18 Xiangtai road, Liuyang Economic Development zone, Changsha, Hunan Province, China  SRN: CN-MF-000042754  Tel :86-731-83167851  Website: <a href="http://www.gmmmedical.net">www.gmmmedical.net</a> Email: <a href="mailto:sales@greatmade.com.cn">sales@greatmade.com.cn</a></p>
	<p>Name : SUNGO Europe B.V.  Add: Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands  SRN: NL-AR-000000247</p>

## 9. Revision information

Rev.	Revision description	Revised date
A.0	First edition	2025/8/12
A.1	Includes English, German language	2025/8/14
A.2	Includes English, German, French, Spanish, Portuguese language	2025/8/15

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